

# ProductCHANGE Announcement

## Serological Pipettes

**NEW, improved packaging and pipette design. Improved accuracy.**


The Globe Scientific serological pipette packaging has been redesigned with these value added features:

- Heavy-duty carton prevents pipette breakage
- Individually wrapped pipettes are sold in cartons that feature color coded labels for easy size identification and dispenser openings for benchtop convenience

In addition each pipette now features:


- Crisp and bold easy-to-read graduation marks
- Reverse graduations on all volumes
- Negative graduations on all volumes
- 1, 2 and 5mL pipettes feature one-piece construction
- 10, 25 and 50mL pipettes feature three-piece construction

**Globe Scientific serological pipettes are manufactured under strict ISO 9001:2008 standards in a class 100,000 cleanroom.** They are non-pyrogenic, non-cytotoxic, RNase-free and DNase-free. The pipettes are equipped with an aerosol barrier plug that prevents both liquid and aerosols from contaminating the pipettor. Printed black graduations prevent ambiguity and guarantee precise reading with less than 2% inaccuracy at full volume. The pipettes are lot certified for accuracy.



	Overall Pipette Length	
	New Version	Previous Version
1mL	275mm	270mm
2mL	275mm	270mm
5mL	342mm	296mm
10mL	345mm	297mm
25mL	345mm	300mm
50mL	345mm	365mm

*Shorter style pipettes coming soon!*



*Reverse and negative graduations on all volumes*



*Crisp and bold easy-to-read graduation marks*

**Full Feature List:**

- Manufactured under strict ISO 9001:2008 standards
- Produced in a class 100,000 cleanroom
- Non-pyrogenic and non-cytotoxic
- RNase-free and DNase-free
- Sterile pipettes are sterilized by gamma irradiation (SAL 10<sup>-6</sup>)
- Aerosol barrier filter plug
- Reverse graduations
- Negative graduations
- Graduation marks are impervious to alcohol wipes
- Calibrated to deliver (TD) – lot certified for accuracy
- Color coded
- Fits all major brands of pipettors
- Individually wrapped pipettes are packaged in convenient color-coded heavy-duty dispenser boxes

**Products**

Item#	Description	Packaging
1700	1mL, sterile serological pipette individually wrapped	1 per pack, 500 per case
1710	1mL, sterile serological pipette bulk packed	25 per bag, 1000 per case
1715	1mL, non-sterile serological pipette bulk packed	25 per bag, 1000 per case
1720	2mL, sterile serological pipette individually wrapped	1 per pack, 500 per case
1730	2mL, sterile serological pipette bulk packed	25 per bag, 1000 per case
1735	2mL, non-sterile serological pipette bulk packed	25 per bag, 1000 per case
1740	5mL, sterile serological pipette individually wrapped	1 per pack, 250 per case
1750	5mL, sterile serological pipette bulk packed	25 per bag, 375 per case
1755	5mL, non-sterile serological pipette bulk packed	25 per bag, 375 per case
1760	10mL, sterile serological pipette individually wrapped	1 per pack, 200 per case
1770	10mL, sterile serological pipette bulk packed	25 per bag, 250 per case
1775	10mL, non-sterile serological pipette bulk packed	25 per bag, 250 per case
1780	25mL, sterile serological pipette individually wrapped	1 per pack, 100 per case
1782	25mL, sterile serological pipette bulk packed	25 per bag, 100 per case
1784	25mL, non-sterile serological pipette bulk packed	25 per bag, 100 per case
1820	50mL, sterile serological pipette individually wrapped	1 per pack, 100 per case
1822	50mL, sterile serological pipette bulk packed	25 per bag, 100 per case



*The individually wrapped pipettes are sold sterile in easy-to-open polyethylene/medical grade paper peel packs.*



*Bulk pipettes are packaged in 25-piece bags and are available sterile or non-sterile. (The 50mL bulk pipettes are available in sterile only.)*



[www.globescientific.com](http://www.globescientific.com)

**Globe Scientific Inc.**  
610 Winters Avenue • Paramus, NJ 07652

**Customer service, sales and technical support:**  
Phone 1 (800) 394-4562 • 1 (201) 599-1400  
Fax 1 (201) 599-1406  
E-mail [mail@globescientific.com](mailto:mail@globescientific.com)