

Blood Dilution Vial, Polystyrene (PS)

Item# 107047



Description

These Globe Scientific polystyrene (PS) cell counter cups come with a polyethylene (PE) snap on lid. The cups are graduated at the 10mL and 20mL marks.

Features

- Excellent clarity
- 35mL capacity
- Graduated at 10mL and 20mL
- Snap cap

Technical Specifications

Material – vial: Polystyrene (PS)

Material – cap: Polyethylene (PE)

Volume: 35mL

Dimensions – vial: 31mm Dia. x 55mm

Dimensions – cap: 38mm Dia. x 6mm

Temperature Range – vial: -10°C – +70°C

Temperature Range – cap: -50°C – +75°C

Sterility: Non-sterile

Autoclavable: No

Products

Item#	Description	Packaging
107047	Blood Cell Counting Vial with Attached Snap Cap	250/bag, 4 bags/case (1000 each)

