

## Wide Bore Transfer Pipet, 7mL

Item# 135040

### Technical Specifications

Material: Low Density Polyethylene (LDPE)  
Length: 129mm (+/- 0.5mm)  
Capacity: 7mL  
Sterility: Non-Sterile  
Autoclavable: No  
Temperature Range: -50°C to +75°C

### Chemical Resistance

- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with halogenated hydrocarbons.

### Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

### Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity

