

Exact Volume Transfer Pipet, 50 μ L (0.05mL)

Item# 139110

Technical Specifications

- Material: Low Density Polyethylene (LDPE)
- Length: 59mm (+/- 1mm)
- Capacity: 50 μ L
- Bulb Draw: 50 μ L
- Sterility: Non-Sterile
- Autoclavable: No
- Temperature Range: -50°C to +75°C

Chemical Resistance

- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with halogenated hydrocarbons.

Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity

