

ALL NEW!

globe[®]
SCIENTIFIC INC.



24 Hour Urine Collection Containers

New 24-hour urine collection containers offer sizes, configurations, and features required to meet most 24-hour collection needs. Upright configurations for efficient storage and low form designs for stability are all designed with easy to read, molded graduations and leak resistant closures. Amber polyethylene construction protects samples from light and molded handles make handling easy. The new TransferTop™ option includes an integrated mechanism for “closed system” transfer to vacuum urine sample tubes.

- Amber polyethylene (PE) construction
- Easy to read molded graduations
- Molded handles for easy handling
- Leak resistant closures
- TransferTop™ option for improved lab safety and hygiene

Item #	Description	Unit
108010	3 Liter, with attached TransferTop™ screwcap, amber	30
108020A	2 Liter, with attached screwcap, amber	54
108025A	2.5 Liter, with attached screwcap, amber	30
108030A	3 Liter, with attached screwcap, amber	30

108010 — 3L TransferTop™ 24 Hour Urine Collection Container features an integrated mechanism for “closed system” transfer of urine directly from the container to vacuum tubes.

globe[®]
SCIENTIFIC INC.

www.globescientific.com

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

Certified
SBA WOSB
Woman Owned Small Business

certified
WBENC
WOMEN'S BUSINESS ENTERPRISE