

# The best GHS compliant wash bottle ever!

**COMPLIES WITH THE LATEST GHS AND OSHA REQUIREMENTS FOR 2015**



**NEW!**

Bottles meet the latest requirements of GHS & OSHA HAZCOM 29 CFR 1910.1200 for workplace labels

*Sold by the each or 5 pack!*

- All bottles include Globally Harmonized Pictograms, NFPA diamond, risk phrases and safety information
- The only GHS wash bottle to be printed in 3 languages — (English, French and Spanish)
- The only GHS wash bottle to come fully assembled
- Safety vented
- Patented DripLok® valve prevents dangerous solvent drips\*
- The specially designed jet-tip provides precise dispensing control
- Wide mouth opening for easy filling
- Color coded

Item #	Capacity	Solvent	Color	Unit(s)
WGW537VTML-GHS	500mL	Acetone	Red	Each/5 pk
WGW538VTML-GHS	500mL	Methanol	Green	Each/5 pk
WGW539VTML-GHS	500mL	Isopropanol	Yellow	Each/5 pk
WGW541VTML-GHS	500mL	Ethanol	Orange	Each/5 pk
WGW601PML-GHS*	500mL	Distilled Water	White	Each/5 pk
WGW500AST-GHS	500mL	Assorted (1 Each)	—	5 pk

\*Distilled water bottles are not safety vented.

**globe**<sup>®</sup>  
SCIENTIFIC INC.

Customer service, sales and technical support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400 • Fax: 1 (201) 599-1406

E-mail: [mail@globescientific.com](mailto:mail@globescientific.com) • [www.globescientific.com](http://www.globescientific.com)

© Copyright 2015 Globe Scientific Inc. All rights reserved. GHSwashbottle 15056


Bottles meet the latest requirements of GHS & OSHA HAZCOM 29 CFR 1910.1200 for workplace labels

Specially designed jet tip provides precise dispensing control

Color coded

**No assembly required!**

Bottles arrive fully assembled and ready for use — by the each or neatly shrink-wrapped in packs of 5



Safety vented: patented DripLok® valve prevents solvent drips

Chemical name in three languages

Chemical formula and CAS (Chemical Abstract Service) reference number

GHS Pictograms and signal words in three languages



Wide mouth for easy filling

Additional safety pictograms and information: required protective equipment and target organs in three languages.

U.S. standard NFPA (National Fire Protection Association) diamond

## Simple and safe transport of up to 6 wash bottles.

### Polypropylene, Cluster Design

- Strong lightweight construction
- Will accept most 500mL bottles
- Internal cavity diameter of 90mm

Item #	Description	Unit
BGS024	Bottle carrier with handle, holds 6 x 500mL bottles	Each

### Polyamide Coated Steel

- Offers better heat and chemical resistance than other coatings

Item #	Description	Unit
BGS040	6-place bottle carrier with handle, for 60mm bottles	Each
BGS042	6-place bottle carrier with handle, for 74mm bottles	Each



Polypropylene bottle carrier



Polyamide coated steel bottle carriers

**globe**<sup>®</sup>  
SCIENTIFIC INC.

Customer service, sales and technical support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400 • Fax: 1 (201) 599-1406

E-mail: mail@globescientific.com • www.globescientific.com