

GTR-HA

Analog Tube Rotator



*GTR-HA Analog Tube Rotator shown
with 15mL and 50mL tube clips*

The GTR-HA Tube Rotator provides slow to vigorous mixing of your biological samples and solutions in a simple, lightweight and easy to use design. Simply turn on the power and adjust the speed using the analog dial for your desired level of mixing. The unit includes a rotisserie-style stainless steel platen and interchangeable tube holder clips for 1.5mL, 15mL and 50mL tubes.

Features:

- End over end motion/agitation
- Provides gentle and effective mixing of biological samples
- Adjustable speed range from 0 to 80rpm
- Unit includes an assortment of twelve 15mL, twelve 50mL and thirty-six 1.5mL/2.0mL tube holder clips
- Simple, robust design for years of worry-free operation

Analog Tube Rotator



Simple analog dial for easy speed adjustment



Tube clips are easily changed by sliding into the channels on the rotisserie tube rack



1.5/2.0mL, 15mL and 50mL tube clips are included and can be mixed and matched on the rotisserie tube holder

GTR-HA Specifications

Tube Capacity	Forty-eight 1.5mL tubes, sixteen 15mL tubes or eight 50mL tubes
Speed Range	0 - 80rpm
Operation Mode	Continuous
Voltage	100-240V, 50/60Hz
Power	20W
Dimensions (WxDxH)	300 x 220 x 310mm (11.8 x 8.7 x 12.2in)
Weight	2.7kg (6.0lbs)
Operating Temperature	5° - 40° C
Permissible Relative Humidity	80% RH (non-condensing)
Protection Class	IP21
Warranty	2 Years

ORDERING INFORMATION

Item #	Description	Unit
GTR-HA	Analog Tube Rotator, Horizontal, Adjustable Speed, 120-240v, 50/60Hz, (Includes 36 each 1.5mL/2.0mL, 12 each 15mL and 12 each 50mL Tube Holder Clips)	Each
GTR-HA-1.5	Tube Holder Clips for 1.5mL/2.0mL Microcentrifuge Tubes	12
GTR-HA-15	Tube Holder Clips for 15mL Centrifuge Tubes	12
GTR-HA-50	Tube Holder Clips for 50mL Centrifuge Tubes	12