

Get a Grip with the GRIP RACK!

"All-in-one" convenience conserves valuable benchtop real estate!

These racks feature flexible rubber grips that firmly hold a variety of tube sizes in a single rack.



Low profile 456920 and 456921 allow for easy viewing of tube labels.

Fully Autoclavable

Innovative Grip Racks

These specially designed test tube racks feature rubber grips that hold test tubes in place when decanting liquid (photo A). Item #456925 & #456926 feature a unique tube ejector that allows the user to easily eject all tubes into a waste receptacle (photos B and C).

- 50-place (5 x 10 configuration)
- Single rack accommodates a variety of tube sizes
- Holds tubes securely while decanting

Low Profile Grip Rack

Holds various size tubes up to 17mm diameter. Dimensions: 254 x 140 x 51mm.

- Low profile rack allows for easy viewing of labels on tubes

Item #	Description	Unit
456920	Autoclavable grip rack, blue	Each
456921	Autoclavable grip rack, orange	Each

Grip Rack with Tube Ejector


Holds various size tubes up to 16mm diameter. Dimensions: 260 x 140 x 82mm.

- Unique tube ejector for simultaneous ejection of all tubes


Item #	Description	Unit
456925	Autoclavable grip rack, blue	Each
456926	Autoclavable grip rack, orange	Each



Both racks feature rubber grips that hold tubes securely in place when decanting liquids.



Push up to eject tubes



Items 456925 & 456926 feature a unique tube ejector for easy disposal of tubes.

globe[®]
SCIENTIFIC INC.

Customer service, sales and technical support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400 •

Fax: 1 (201) 599-1406 • E-mail: mail@globescientific.com

© Copyright 2013 Globe Scientific Inc. All rights reserved. grprk 13160