

Ice Buckets & Trays

Globe Scientific ice buckets and ice trays are manufactured from polyurethane foam for excellent insulating properties and good chemical resistance.

Rectangular Ice Buckets

Globe Scientific rectangular ice buckets are suitable for use with liquid nitrogen, dry ice, wet ice, and salt solutions. Lids included with ice buckets and select ice trays limit evaporation for extended chilling. Molded handles provide a secure grip for easy transport.



Rectangular Ice Buckets

Volume	Dimensions	Unit	Item#			
			Black	Blue	Green	Red
2.5 Liter	340mm x 290mm x 120mm	Each	455010K	455010B	455010G	455010R
4.5 Liter	330mm x 280mm x 185mm	Each	455015K	455015B	455015G	455015R

All NEW!

Ice Buckets & Trays

Globe Scientific ice buckets and ice trays are manufactured from polyurethane foam for excellent insulating properties and good chemical resistance.

Rectangular Ice Trays

Globe Scientific rectangular ice trays are suitable for use with liquid nitrogen, dry ice, wet ice, and salt solutions. Lids included with select ice trays limit evaporation for extended chilling. Molded handles provide a secure grip for easy transport.



Ice Trays Without Lid

Volume	Dimensions	Unit	Item#			
			Black	Blue	Green	Red
1 Liter	210mm x 190mm x 105mm	Each	455020K	455020B	455020G	455020R
4 Liter	380mm x 270mm x 152mm	Each	455022K	455022B	455022G	455022R
9 Liter	500mm x 360mm x 160mm	Each	455024K	455024B	455024G	455024R

Ice Trays With Lid

Volume	Dimensions	Unit	Item#			
			Black	Blue	Green	Red
1 Liter	210mm x 190mm x 105mm	Each	455021K	455021B	455021G	455021R
4 Liter	380mm x 270mm x 152mm	Each	455023K	455023B	455023G	455023R
9 Liter	500mm x 360mm x 160mm	Each	455025K	455025B	455025G	455025R



www.globescientific.com

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

