

NEW!

globe[®]
SCIENTIFIC INC.

Diamond High-Performance Pipettor Starter Packs



EXTRAORDINARY VALUES!

*All the features and benefits of Diamond Advance™ Pipettors, Diamond Pro™ Pipettors, and Globe Scientific Certified Pipette Tips **PLUS:***

Efficiency –

Each starter pack includes four sizes of pipettors enabling volumes from 0.5uL to 1000uL, a pipet stand, and a **full box of tips for each size!** That's 36 racks of tips... enough to really get work done.

Savings –

Four high performance pipettors, the stand, and tips for less than the normal price of the pipettors alone!

Choice –

Starter packs are available with either Diamond Advance™ Ultralight High Performance Pipettors or Diamond Pro™ High Performance Pipettors to meet nearly all performance needs and price points.

Quality –

Starter packs include Globe Scientific Certified pipet tips that are certified sterile, RNase, DNase, Pyrogen, DNA, and PCR inhibitor free making them suitable for the most critical research.

Accuracy –

Pipettors are calibrated in an ISO 17025 accredited laboratory and both Diamond Advance™ and Diamond Pro™ pipettors have padded handles to reduce heat transfer for consistent accuracy.

DiamondAdvance™

ultralight high-performance pipettor **STARTER PACK**

over
45% OFF
list price of items if
purchased separately!

Each Diamond Advance Starter Pack contains:

- **4 Adjustable volume pipettors:**
 - 0.5-10uL Diamond Advance Pipettor (1)
 - 2-20uL Diamond Advance Pipettor (1)
 - 20-200uL Diamond Advance Pipettor (1)
 - 100-1000uL Diamond Advance Pipettor (1)
- **4 Full boxes of racked, sterile, certified*, low-retention pipette tips:**
 - 0.1-10uL (1 Box of 960 tips, 10 racks per box)
 - 0.1-20uL (1 Box of 960 tips, 10 racks per box)
 - 1-200uL (1 Box of 960 tips, 10 racks per box)
 - 100-1250uL (1 Box of 576 tips, 6 racks per box)
- **6-Place carousel pipette stand**



Item #3341-COMBO

**Certified sterile, RNase, DNase, Pyrogen, DNA, and PCR inhibitor free.*

DiamondPRO™

high-performance pipettor **STARTER PACK**

over
45% OFF
list price of items if
purchased separately!

Each Diamond Pro Starter Pack contains:

- **4 adjustable volume pipettors:**
 - 0.5-10uL Diamond Pro Pipettor (1)
 - 2-20uL Diamond Pro Pipettor (1)
 - 20-200uL Diamond Pro Pipettor (1)
 - 100-1000uL Diamond Pro Pipettor (1)
- **4 full boxes of racked, sterile, certified*, general purpose pipette tips:**
 - 0.1-10uL (1 Box of 960 tips, 10 racks per box)
 - 0.1-20uL (1 Box of 960 tips, 10 racks per box)
 - 1-200uL (1 Box of 960 tips, 10 racks per box)
 - 100-1250uL (1 Box of 576 tips, 6 racks per box)
- **4-place pipette stand**



Item #3311-COMBO

**Certified sterile, RNase, DNase, Pyrogen, DNA, and PCR inhibitor free.*

globe®
SCIENTIFIC INC.

Customer service, sales and technical support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400 • Fax: 1 (201) 599-1406

E-mail: mail@globescientific.com • www.globescientific.com

© Copyright 2016 Globe Scientific Inc. All rights reserved. pkits 16354