


Straight and Curved **Transfer Syringes**

Globe Scientific's 12cc (12mL) Transfer Syringes are offered in both straight tip and curved tip for multiple uses and methods. The syringes are made from medical grade polypropylene. Printed graduations are bold and easy to read. The plunger can be removed for cleaning. Typical uses include oral irrigation, animal feeding and more.

Features:

- Manufactured from medical grade polypropylene (PP)
- Bold, printed graduations to 12mL at 1mL increments
- 12cc (12mL) volume working capacity
- Removable plunger for easy cleaning



NEW!

Item #	Description	Unit
159012	Transfer Syringe, 12cc (12mL), Printed Graduations, Straight Tip, 50/bag	50
159014	Transfer Syringe, 12cc (12mL), Printed Graduations, Curved Tip, 50/bag	50