

1.5mL Sample Cup, Polystyrene (PS)

Item# 110610



Description

These Globe Scientific multi-purpose sample cups are designed for use on the most popular instruments in the industry. The sample cups are made to exact specifications and provide excellent optical properties.

Features

- 1.5mL capacity
- Excellent clarity

Technical Specifications

Material: Polystyrene (PS)

Dimensions: 13.8mm x 22.60mm

Autoclavable: No

Brittleness temperature: 20°C

Maximum temperature: 90°C

Sterility: Non-sterile

Products

Item#	Description	Unit
110610	1.5mL sample cup, PS	1000

