

0.25mL Sample Cup, Polystyrene (PS)

Item# 110711



Description

These Globe Scientific multi-purpose sample cups are designed for use on the most popular instruments in the industry. The sample cups are made to exact specifications and provide excellent optical properties.

Features

- 0.25mL capacity
- Excellent clarity

Technical Specifications

Material: Polystyrene (PS)

Dimensions: 14mm x 16mm

Autoclavable: No

Temperature range: -10°C to +70°C

Sterility: Non-sterile

Products

Item#	Description	Unit
110711	0.25mL sample cup, PS	1000

