

16 x 100mm (12mL) Polypropylene (PP) Test Tube

Item# 111040



Description

These Globe Scientific test tubes are molded from premium grade polypropylene (PP) for impact strength and stability to common acids, alkalies, and solvents. The tubes have a maximum spin rate of 3000 x g and are autoclavable and freezable.

Features

- 12mL capacity
- 16 x 100mm
- Graduated at 2.5mL, 5mL & 10mL
- With rim

Technical Specifications

Material: Polypropylene (PP)
 Spin Rate: 3000 x g
 Autoclavable: Yes
 Brittleness temperature: 0°C
 Maximum Temperature: 135°C
 Sterility: Non-sterile

Products

Item#	Description	Unit
111040	16 x 100mm, 12mL test tube, PP	2000

