


12mL Conical Bottom Centrifuge Tube

Item# 112010

Description

These Globe Scientific general purpose centrifuge tubes feature conical bottoms and are made from clear polystyrene (PS). The tubes are graduated at 0.5, 1, 2.5, 5 and 10mL.



Features

- 12mL capacity
- 16 x 100mm
- Conical bottom
- Polystyrene (PS)
- Graduated to 10mL
- Graduated at 0.5mL, 1mL, 2.5mL, 5mL and 10mL

Technical Specifications

Spin Rate: 1500 RCF
 Autoclavable: No
 Brittleness temperature: 20°C
 Maximum temperature: 90°C
 Sterility: Non-sterile

Products

| Item# | Description | Units |
|--------|------------------------------------------------|------------|
| 112010 | 12mL centrifuge tube, PS, 500/bag, 4 bags/case | 500 & 2000 |