

1-200µL Certified Pipette Tips

Item# 151150, 151150R-96, 151150RS-96, 151150RF, 151150Y, 151150YR-96, 151150YRS-96, 151150YRF



151150



151150Y

Description

Globe Scientific CERTIFIED pipette tips are certified free from Human DNA, DNase, RNase, Pyrogens and PCR Inhibitors. These premium quality pipette tips are manufactured in a cleanroom using fully automated equipment. They are produced from medical grade virgin polypropylene (PP). All are graduated for accurate measurements. The universal fit design provides an airtight fit on most brands of pipettors. They are available in unique self-standing bags, racks and reloading stacks.

Features

- Produced on state-of-the-art, fully automated production equipment
- Certified free from Human DNA, DNase, RNase, Pyrogens and PCR Inhibitors
- Heavy-metal free
- USP Class IV and ISO 10993 certified
- Available non-sterile in bags, sterile or non-sterile in racks, and non-sterile in reloading stacks

Packaging options:



Bulk in bag

Racked

Reloading stacks

Products


Item #	Volume Range	Tip Length	Tip Color	Sterility	Packaging	Unit
151150	1 - 200µl	54mm	Natural	–	Bulk in Bag	1000 Tips/Bag
151150R-96	1 - 200µl	54mm	Natural	–	Racked	96 Tips/Rack, 10 Racks/Box
151150RF	1 - 200µl	54mm	Natural	–	Reloading Stack	96 Tips/Tray, 10 Trays/Stack
151150RS-96	1 - 200µl	54mm	Natural	Sterile	Racked	96 Tips/Rack, 10 Racks/Box
151150Y	1 - 200µl	54mm	Yellow	–	Bulk in Bag	1000 Tips/Bag
151150YR-96	1 - 200µl	54mm	Yellow	–	Racked	96 Tips/Rack, 10 Racks/Box
151150YRF	1 - 200µl	54mm	Yellow	–	Reloading Stack	96 Tips/Tray, 10 Trays/Stack
151150YRS-96	1 - 200µl	54mm	Yellow	Sterile	Racked	96 Tips/Rack, 10 Racks/Box

1-200 μ L Certified Pipette Tips

Item# 151150, 151150R-96, 151150RS-96, 151150RF, 151150Y, 151150YR-96, 151150YRS-96, 151150YRF

Technical Specifications

Material: Virgin Polypropylene (PP)
Autoclavable: Yes
Graduations: Molded, 10 μ L, 50 μ L, 100 μ L, 200 μ L,
Maximum Volume: 200 μ L
Sterility: Sterile:
151150RS-96, 151150YRS-96
Non-Sterile:
151150, 151150R-96, 151150RF, 151150Y,
151150YR-96, 151150YRS-96



Tips shown actual size



Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400
Fax 1 (201) 599-1406
E-mail mail@globescientific.com

Globe Scientific Inc.

400 Corporate Drive
Mahwah, NJ 07430
www.globescientific.com

