

INTRODUCING!



96 Deep-Well Plate

HTS friendly – compatible with most liquid-handling robots.

Globe Scientific is excited to introduce these new **STERILE** 96 deep-well plates! They are compatible with numerous applications including compound storage, lysing methods, molecular and cell biology and many more. Each round well incorporates highly polished walls, and a U-shaped bottom for maximum sample recovery and centrifugation. The SBS footprint allows use with automated liquid handling robots. These deep-well plates are manufactured from virgin polypropylene to ensure purity and sample integrity, chemical resistance and compatibility with long-term storage needs. They are available in cases of 50 plates (10 packs of 5 plates)



Features:

- Certified RNase, DNase free
- Sterile
- SBS footprint — meets ANSI/SLAS 1-2004
- Flatness tested for reliable performance on automated equipment
- Smooth well surfaces ensure maximum sample recovery
- Centrifugation to 4,000 x g
- Raised rim facilitates secure sealing with films or sealing mats
- High clarity for easy sample viewing
- Easy-to-read alpha-numeric indexing
- Notched corners for easy orientation
- Autoclavable
- Stackable
- Recyclable

Item #	Description	Unit
V-DWP-96-U-S	Deep Well Plate, 96 Well, Round Wells, U-Bottom, STERILE, 2.0mL, 5 Plates/Pack, 10 Packs/Case	50

Outstanding performance with top-notch features and benefits . . . choose Globe's new 96 well deep well plate!



www.globescientific.com

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com